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The screenshot displays the GCMS Real Time Analysis software interface. The title bar reads "GCMS Real Time Analysis (Admin) - [Acquisition - Disabled, 060918-3.qgd(Live1), 060909.qgd]". The menu bar includes File, Edit, View, Method, Instrument, Acquisition, Data, Tools, Window, and Help. The File menu is open, showing options such as "New Method File", "Open Method File...", "Close Method File", "Save Method File", "Save Method File As...", "Save Method As Template...", "Load Method Parameters...", "Open Reference Data File...", "Close Reference Data File", "Select Project(Folder)...", "File Search...", "System Administration...", "Audit Trail...", "Print Setup...", "Print Method File", "Method File Properties...", and a list of recent method files. The main window shows a chromatogram with a yellow background and a grid. The x-axis is labeled "Time: 0:00 / 10:00 min" and "Scan #: 0". The y-axis is labeled "Intensity: 0" and "Group: 1". The baseline is flat at "Base Peak Intensity: 0". The x-axis has major ticks at 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, and 8.5. The y-axis has a major tick at 50. The bottom right corner shows a small inset chromatogram with a y-axis labeled "Split" and a x-axis labeled "min".

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The screenshot displays the GCMS Real Time Analysis software interface. The top menu bar includes File, Edit, View, Method, Instrument, Acquisition, Data, Tools, Window, and Help. The toolbar contains various icons for file operations and analysis. The main window is divided into several sections:

- Acquisition Panel (Left):** Contains buttons for Top, Method, Sample Logon, Standby, Start, Stop, Method Detail, Snapshot, and Data.
- Chromatogram (Top Center):** Shows a plot of intensity versus time (min) with a yellow background. The x-axis ranges from 0.0 to 45.0 minutes. The y-axis ranges from 0.0 to 2.0. A single peak is visible at approximately 15.0 minutes.
- Parameters Panel (Bottom Left):** Displays various acquisition parameters:
 - Inj. Port: SP12-PT
 - Inj. Mast Port: (blank)
 - Column Oven Temp: 50.0 °C
 - Injection Temp: 200.0 °C
 - Injection Mode: Split
 - Split Ratio: 10.0
 - Flow Control Mode: Linear Velocity
 - Pressure: 300.2 kPa
 - Total Flow: 25.0 mL/min
 - Column Flow: 3.61 mL/min
 - Linear Velocity: 40.3 m/sec
 - Sample Flow: 3.0 mL/min
- Program Panel (Bottom Right):** Shows the temperature program settings:
 - Program: Column Oven Temperature
 - Table with columns: Run, Final Temperature, Hold
 - Total Program Time: 40.37 min

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File Edit View Method Instrument Acquisition Data Tools Window Help

Acquisition

Sample Name : 080918-3
Sample ID :
Data Description :

Line1-MS Rel. Time 0:00:49.00 min Scan 0-1 Intensity 0 Group 1 C: 9.77

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